



PROTECTIVE COVERALLS TECHNICAL DATE SHEET

Protective Coverall designed to help protect against liquid chemicals (Type 3 - Type 4) hazardous dust (Type 5), certain light liquid splashes (Type 6)

Three panel hood design on coveralls for compatibility with complementary PPE

Seamless shoulder and sleeve tops for fewer entry points for contaminants and increased comfort

Elasticized waist and ankles and knit or elastic cuffs on all models for increased comfort and mobility

Two-way zipper for on/off convenience



Model: OLEY0909

Model: OLEY2311

Sewn from Laminated Non-Woven Fabric with High Quality, Light Weight and Easy to Wear. Made of Soft, Sturdy, Breathable and Waterproof Fabric.

To Increase Protection Against Less Dangerous Chemicals and Dusts; Two-layer, Hooded and Covered Zipper is used with Laminated Fabric.

Made from 57 gsm laminated material.

EN ISO 13688 Protective Clothing - General Requirements

EN 13034+ A1 Protective Clothing Against Liquid Chemicals (Type 6-B)

EN ISO 13982-1 Protective Clothing Against Solid Particules (Type 5-B)


EN 14605+A1 Protective Clothing Against Liquid Chemicals (Type 3-B, Type 4-B)

EN 14126 Protective Clothing Performance Requirements and Test Methods For Protective Clothing Against Infective Agents




Liquid Protection

Type 6 - Light Liquid Splashes
(EN 13034:2005 + A1:2009)




Liquid Protection

Type 4 - Liquid Spray
(EN 14605:2005 + A1:2009)




Dust Protection

Type 5 - Solid Airborne Particulates
EN ISO 13982-1:2004




Anti-static

Anti-static coating* on both sides
(EN 1149-5:2008)



Nuclear

Radioactive particulates Class 2
(EN 1073-2:2002 except EN863 puncture)




Biohazard

Types 4-B, 5-B & 6-B - Infective Agents
(EN 14126:2003 except ISO 16604)



* All apparel must be suitably grounded for anti-static treatment to be effective.
Does not offer protection against radiation.



Breathable Fabric PP+PE



Limited use



Do not dry



Do not wash



Do not dry-clean



Do not bleach




Do not iron




Flammable



Superior Protection




Verstärkte Schrittbereichstützung
Reinforced Crotch Support



Personal Protective Coverall
Persönlicher Schutzoverall

- S
- M
- L
- XL
- 2XL
- 3XL
- 4XL
- 5XL



1 STCK.
PCS.



COVERALLS SIZE CHART

Size	Chest	Length
M-L	92-108 cm	170-180 cm
L-XL	108-116 cm	180-190 cm
XL-XXL	116-124 cm	190-195 cm
3XL-4XL	124-135 cm	190-195 cm


Label Example

PROTECTIVE COVERALLS
OLEY0909 / OLEY2311




8 6 8 2 4 5 0 8 2 0 3 9 5

Production Date: :11/2020
Expiry Date: :11/2025
LOT No : :202101/102
Size : :2XL



Manufacturer: Oley Tıbbi Ürünler Medikal Mlz. San. Tic. Ltd Şti

Address: Parseller Mh. Menekşe Cd. Feyyaz Sk. No:2/3
Ümraniye İstanbul TURKEY www.oleytibbiurunler.com
info@oleytibbiurunler.com



Packing List

1 Master Box can carry 50 pcs.(From L size to 3XL Size)
Box Dimensions are 60*40*40 cm
1 box is 9,32 kg

-Packing List may change based on sizes, please indicate the sizes