



# OLEY NON STERILE STANDARD GOWNS TECHNICAL DATA SHEET



Designed specifically for non-sterile procedures, this gown is constructed from SMS fabric and features a tie neck closure, knit cuffs, and a waist belt.

The materials used in the gowns is made of medical grade materials.


Gowns are in closed packages one by one.

The knit cuffs of the gowns are thin and tight for ease of use.

35-40 gsm fabric thicknesses

CE Marked

# OLEY NON STERILE GOWNS



Production Date: :10/2020  
Expiry Date: :10/2023  
LOT No :202009/368  
Size : L  
40 gsm

Manufacturer: Oley Tıbbi Ürünler Medikal  
Mlz. San. Tic. Ltd Şti  
Address: Parseller Mh. Menekşe Cd. Feyyaz Sk.  
No:2/3 Ümraniye İstanbul TURKEY  
www.oleytibbiurunler.com  
info@oleytibbiurunler.com



Gowns are packed one by one in polybags. 1 master box contains 60 pcs.  
Dimensions of Master Box:60\*40\*40 cm



Limited use



Do not dry



Do not wash



Do not dry-clean



Do not bleach



Do not iron



Flammable

## Accessory Storage and Disposal

- Store in dry, clean conditions in original packaging, away from direct sunlight, sources of high temperature, and solvent vapors
- Store within the temperature range -20°C to +25°C (-4°F to +68°F) and with relative humidity below 80%
- Expected shelf life is three years from date of manufacture when stored as stated
- Replace product if damaged, heavily contaminated or in accordance with local regulations.
- Handle and dispose of contaminated product with care and in accordance with applicable regulations